

# Vent II Introduction to Ventilator Management

- Medical Practice Improvement Project
- Paul Wisniewski, DO
- January 9, 2025



# Disclosures

- None



# Learning Objectives

- 1. Learning the difference between time cycle and volume cycled ventilation
- 2. How to adjust the ventilator to compensate for high PCO<sub>2</sub>
- 3. Management strategies improve oxygenation
- 4. Evaluation of ABG



# Modes of Ventilation

- Volume-cycled ventilators:
  - Gas flows to the patient until a preset volume is delivered to the ventilator circuit, even if this entails a very high airway pressure.
  - Individuals require mechanical ventilation for different reasons.
- Examples:
  - Assist Control ( A/C)
  - Synchronized intermittent mandatory ventilation ( SIMV)



# Modes of Ventilation

- Pressure Control Ventilation (PCV) is a mode of mechanical ventilation where:
- The ventilator controls all phases of the breath based on set parameters
- Breath delivery is at a set inspiratory pressure, optimizing oxygenation and reducing barotrauma risk
- The pressure remains constant, but the volume of air delivered can change based on lung compliance and airway resistance
- It is used to regulate pressures applied during mechanical ventilation



# How to Correct CO<sub>2</sub> and PaO<sub>2</sub>

## Oxygenation

- Increased by increasing FiO<sub>2</sub>
- Mean Airway Pressure
- Peep

## Carbon Dioxide Clearance

- Increasing TV
- Increasing Rate
- ( increased Minute Ventilation)



# Setting the Ventilator Settings

- Volume Cycled Ventilation
  - **Tidal Volume** is set at less than 10 ml/kg for ideal body wt
    - Lung protective strategy is 6-8 ml/kg for ideal body wt
  - **Respiratory rate** : 10-20 BPM
  - **Fio2** starts at 100% and titrates down
  - **Peep** Starts at 5 cm of water and titrated up
  - **Pressure support** starts at 10 cm of water and titrates up for SIMV Ventilation



# SIMV VS A/C = Same mode of ventilation at High Rate

- Example of settings for 75 Kg person: A/C
- R: 14
- Tidal volume 550ml
- Fio2 100%
- PEEP 5 CM H2O

- Example of settings for 75 Kg person: SIMV
- R: 14
- Tidal volume 550ml
- Fio2 100%
- PEEP 5 CM H2O
- Pressure Support 10 cm H2O





# A/C and SIMV Differ when it comes to Spontaneous Breathing

- Patient can breathe spontaneous on either mode
- Rate set at 10 BPM they can breathe at 20 BPM
- The first 10 BMP are the same the difference is in the spontaneous breathes.
- A/C every breathe is the same
- SIMV the spontaneous breathes are based on what the patient can do with pressure support help. PS = IPAP ( BIPAP)



# Adjusting PS to Obtain Desired Tidal Volume



# A/C Drawbacks

- If the patient is Tachypneic they can blow down their CO<sub>2</sub> and become alkalotic.
- Why?
  - High minute ventilation
  - Every Breathe delivers the entire tidal volume
  - Min vent = RR \*TV



# SIMV Drawback

- If the PS is not titrated up appropriately the patient will have
  - High PCOs from low tidal volume and decreased minute ventilation
  - Acidosis from hypoventilation
  - They will attempt to compensate by increasing the rate to blow off CO<sub>2</sub>
- Why?
  - Low minute ventilation
  - Spontaneous Breaths with deliver only what the patient can draw in
  - $\text{Min vent} = \text{RR} * \text{TV}$  ( If they can not get volume they will try to make if up with rate)
  - People that fail SIMV from tachypnea have not had adequate PS added
  - PS does not add to the mean airway pressure. You can titrate as high as needed to get desired tidal volume. ( 10-25 cm H<sub>2</sub>O)



# Carbon Dioxide Calculation

## Adjusting PaCO<sub>2</sub> Using Tidal Volume

$$VT_{new} = \frac{VT_{current} \times PaCO_{2\,current}}{PaCO_{2\,desired}}$$

### Example

- VT current = 400 mL
- PaCO<sub>2</sub> current = 60 mmHg
- PaCO<sub>2</sub> desired = 40 mmHg

$$VT_{new} = \frac{400 \times 60}{40} = 600 \text{ mL}$$

## Adjusting PaCO<sub>2</sub> Using Respiratory Rate

$$RR_{new} = \frac{RR_{current} \times PaCO_{2\,current}}{PaCO_{2\,desired}}$$

### Example (Rate Goes Down)

- RR current = 20
- PaCO<sub>2</sub> current = 30 mmHg
- PaCO<sub>2</sub> desired = 40 mmHg

$$RR_{new} = \frac{20 \times 30}{40} = 15$$



# Adjusting to compensate for high PCO<sub>2</sub>

- Min vent = RR \*TV
- Adjust based on Volume
- Adjust based on Rate
- Example
- PCO<sub>2</sub> is 55mm Hg
- TV 350 ml
- RR 10 BPM
- Desired 40 mmHg PCO<sub>2</sub>



# Calculations

- Tidal Volume
- Equation
- Tidal Volume PCO2 =

$$(\text{PCO2}_{(\text{Have})}) * (\text{TV})_{(\text{Have})}$$

---

$$(\text{PCO2})_{(\text{Want})}$$

- Respiratory Rate equation
- Equation
- RR PCO2 =

$$(\text{PCO2}_{(\text{Have})}) * (\text{RR})_{(\text{Have})}$$

---

$$(\text{PCO2})_{(\text{Want})}$$



# Calculations

- Tidal Volume
- Equation
- Tidal Volume PCO2 =

$$(55) * (350)$$

---

$$(40)_{(Want)}$$

$$= 481.25 = 482 \text{ml}$$

- Respiratory Rate equation
- Equation
- RR PCO2 =

$$(55) * (10)_{(Have)}$$

---

$$(40)_{(Want)}$$

$$= 13.75 = 14 \text{ BPM}$$



# Carbon Dioxide Calculation


## Clinical Pearl

- Increasing PaCO<sub>2</sub> → **decrease VT or RR**
- Decreasing PaCO<sub>2</sub> → **increase VT or RR**
- Assumes stable dead space and CO<sub>2</sub> production



# Management Strategies Improve Oxygenation


- Recruit lung to save the day
- Increased  $FiO_2$  is the rescue maneuver
- If you increase  $FiO_2$  you should increase PEEP
- Need to fix the Right to Left Shunt




# Management Strategies Improve Oxygenation

## Oxygenation

- Increased by increasing  $F_iO_2$
- Mean Airway Pressure
- PEEP



Atelectasis No gas exchange



Interstitial edema take longer for the gas to diffuse over the distance  
It is to the 4<sup>th</sup> power of the distance



# Management Strategies Improve Oxygenation

**High & Low PEEP tables from ARDSnet**

FiO <sub>2</sub>	Low PEEP	High PEEP
0.3	5	5-14
0.4	5-8	14-16
0.5	8-10	16-20
0.6	10	20
0.7	10-14	20
0.8	14	20-22
0.9	14-18	22
1.0	18-24	22-24

PEEP tables don't need to be followed precisely, but can be useful as a general guide. The WHO recommends using a high-PEEP strategy, which seems consistent with available experience thus far with COVID-19. If high PEEPs are used, make sure to keep tidal volumes low to prevent excessively high plateau pressures. APRV is an alternative strategy which would likewise provide high mean airway pressures.

-The Internet Book of Critical Care, @PulmCrit

## PEEP

- As FiO<sub>2</sub> increases, PEEP should also increase

**Lower PEEP/higher FiO<sub>2</sub>**

FiO <sub>2</sub>	0.3	0.4	0.4	0.5	0.5	0.6	0.7	0.7
PEEP	5	5	8	8	10	10	10	12

FiO <sub>2</sub>	0.7	0.8	0.9	0.9	0.9	1.0
PEEP	14	14	14	16	18	18-24

ARDSnet. NEJM 2004; 351, 327



# Management Strategies Improve Oxygenation

<b>FiO2 %</b>	<b>30</b>	<b>40</b>	<b>50</b>	<b>60</b>	<b>70</b>	<b>80</b>	<b>90</b>	<b>100</b>
PEEP (cm H2O)	5	5	8	10	12	14	16	18

**Wisniewski Simplified Version**



# Mean Airway Pressure

- Think of PEEP in terms of the Diastolic BP
- PEEP is the DP
- Mean Arterial Pressure is defined as
- Biggest Contributor to MAP is the DP
- Biggest Contributor to Mean Airway Pressure is PEEP

$$mAP = \frac{1}{3} \text{ systolic pressure} + \frac{2}{3} \text{ diastolic pressure}$$



# Mean Airway Pressure

## Mean Airway Pressure ( $\bar{A}P_{aw}$ ) Equation

$$\bar{P}_{aw} = \left( \frac{T_I}{T_{tot}} \right) P_{IP} + \left( \frac{T_E}{T_{tot}} \right) PEEP$$

### Legend

- $\bar{P}_{aw}$  = Mean airway pressure (cmH<sub>2</sub>O)
- $P_{IP}$  = Peak inspiratory pressure (cmH<sub>2</sub>O)
- $PEEP$  = Positive end-expiratory pressure (cmH<sub>2</sub>O)
- $T_I$  = Inspiratory time (seconds)
- $T_E$  = Expiratory time (seconds)
- $T_{tot}$  = Total respiratory cycle time

$$T_{tot} = T_I + T_E$$



# Mean Airway Pressure

Mean Airway Pressure (A/C Volume)

$$\bar{P}_{aw} = \left( \frac{T_I}{T_{tot}} \right) P_{IP} + \left( \frac{T_E}{T_{tot}} \right) PEEP$$

Common ICU Shortcut (I:E = 1:2)

$$\bar{P}_{aw} = \frac{1}{3} P_{IP} + \frac{2}{3} PEEP$$



# Mean Airway Pressure Increases

## To Increase Mean Pressure

1. increase flow
2. increase PIP
3. increase I:E
4. increase PEEP

- Increasing the Time in Inhalation increase Airway pressure. More time at PIP
- More Flow through the same size tube has to increase pressure
- Increasing the Peak Inspiratory Pressure (PIP) adds to Airway pressure
- PEEP Moves the Lower end of the pressure reading upward



# ABG Analysis

- pH 7.35-7.45
- $\text{PCO}_2$  35- 45 mmHg
- $\text{PO}_2$  60-90 mmHg
- Be +2 to -2
- $\text{HCO}_3$  22-28

We have addressed  
Fixing  $\text{PCO}_2$  ,  $\text{PO}_2$



# ABG Analysis

- Fixing pH
- Acidosis
  - Metabolic
    - BE - < -4
    - pH < 7.35
    - PCO<sub>2</sub> < 35 ( compensating) tachypneic patients are worrisome)
    - HCO<sub>3</sub> < 22
  - Addressed this with Shock lectures  
Resuscitation is the answer
- Respiratory
  - BE- > -4
  - pH < 7.35
  - PCO<sub>2</sub> >45
  - HCO<sub>3</sub> >28 ( Compensatory mechanism)
- Mechanical Ventilation is the answer



# ABG Analysis

- Fixing pH
- Alkalosis
  - Metabolic
    - BE - > +2
    - pH < 7.45
    - PCO<sub>2</sub> >45 ( compensating) slowed breathing (Bradypnea) patients are worrisome)
    - HCO<sub>3</sub> >28
  - Respiratory
    - BE- = neutral
    - pH > 7.45
    - PCO<sub>2</sub> <35
    - HCO<sub>3</sub> = 22 to 28 (slow Compensatory mechanism)



# Alkalosis is Lethal

- A direct relationship between mortality and blood pH exists when the blood pH is  $>7.48$ .
- Mortality rates of 45% and 80% have been noted at blood pH levels of 7.55 and 7.65, respectively .



# Bicarbonate Deficit Calculation

Symbol	Meaning
$\text{HCO}_3^-$ dose	Total sodium bicarbonate needed (mEq)
0.3	Approximate bicarbonate distribution (30% of body weight)
Weight (kg)	Patient body weight in kilograms
Base Deficit	Absolute value of base deficit (mEq/L)

$$\text{HCO}_3^- \text{ dose (mEq)} = 0.3 \times \text{Weight (kg)} \times \text{Base Deficit (mEq/L)}$$

## Worked Example

- Patient: 90 kg
- Base deficit: -10 mEq/L

$$\text{HCO}_3^- = 0.3 \times 90 \times 10 = 270 \text{ mEq}$$

Ampules needed (50 mEq/ampule):

$$270/50 = 5.4 \approx 5-6 \text{ ampules}$$

⚠ Clinically, usually give **half initially**, then reassess ABG.



# Next week

- Look at pressure control ventilation
  - Calculating tidal volume
  - Managing Pressures
  - Increasing oxygenation



# Works Cited

- [Modes Of Mechanical Ventilation Explained | Tracheostomy Education](#)
- [pressure control ventilation – Search](#)
- [simv ventilation with too little ps pressure PV loop - Search Images](#)
- [Picture of Right to left pulmonary shunt - Search Images](#)
- [table for peep and fio2 titration - Search Images](#)
- [table for peep and fio2 titration - Search Images](#)
- [mean arterial pressure equation - Search Images](#)



# Disclaimer for Presentations

The **Presentations** provided by **Cutting Edge Surgical Medical Group**, a division of **Paul J. Wisniewski, DO, Inc.**, are intended to offer general information and guidance based on current research, clinical best practices, and expert opinions in the medical field. These presentations are designed to assist healthcare professionals in making informed decisions regarding patient care, but they are not a substitute for personalized medical advice, diagnosis, or treatment.

## Important Notes:

- The presentations are for informational purposes only and are not intended to replace professional medical judgment. They should be used as a reference and adapted to the specific needs of individual patients.
- Application of these presentations should be made by healthcare providers, taking into account the unique medical history, condition, and circumstances of each patient.
- While **Cutting Edge Surgical Medical Group** strives to provide the most accurate, up-to-date, and evidence-based information, we cannot guarantee that all content on the website is free from errors, omissions, or outdated information. Medical knowledge evolves rapidly, and presentations may be updated periodically.
- **Cutting Edge Surgical Medical Group** does not assume responsibility for the outcomes of any medical decision or intervention based on the use of these presentations. The use of this information is at the user's own discretion.
- Healthcare providers are encouraged to consult the latest peer-reviewed research, professional standards, and individual patient assessments before making clinical decisions.

For specific medical concerns, treatment advice, or patient management, please consult directly with a qualified healthcare provider.

